



| | |
|---------------|----------------|
| Product Code | 26607 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

| | |
|------------|---------------|
| RDM (NSIS) | 2740082 |
| CND | V0402 |
| EAN/UPC | 8023279266078 |
| PARAF | 903162547 |
| GMDN | 12143 |

Description

Instrument tray with angled edges

- With cover
- Dimension: 306 x 196 x h 50 mm
- Autoclavable 121°C

Stainless Steel AISI 304 holloware with minimum thickness 0.6 mm