



Product Code	26610
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2527571
CND	V0402
EAN/UPC	8023279266108
PARAF	903162562
GMDN	42893

Description

S/S AISI 304 - 18/8 KIDNEY DISHES - THICKNESS 0.6 MM - AUTOCLAVABLE 121°C

- SHALLOW KIDNEY DISHES
- Tray 8"- 207x128xh 33 mm 480 ml

Technical Specifications

- Autoclavable 121°C
- AISI 304