

S/S KIDNEY DISH - 280x141x40 mm



Product Code	26612
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2739353
CND	V0402
EAN/UPC	8023279266122
PARAF	903162624
GTIN	08023279266122
GMDN	42893

Description

S/S AISI 304 - 18/8 KIDNEY DISHES - THICKNESS 0.6 MM - AUTOCLAVABLE 121°C

- SHALLOW KIDNEY DISHES
- Tray 12"- 280x141xh 40 mm 1600 ml

Technical Specifications

- Autoclavable 121°C
- AISI 304