

# S/S GALLIPOT WITH LIP diam. 128 mm



Product Code	26618
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2740318
CND	V0402
EAN/UPC	8023279266184
PARAF	903162701
GTIN	08023279266184
GMDN	42893

## Description

S/S AISI 304 - 18/8 GALLIPOTS - AUTOCLAVABLE 121°C

- GALLIPOT WITH LIP
- Ø 128xh 60 mm
- Capacity 450 ml

For all wash Basins and Gallipots we indicate upper diameter