

S/S MEDICINE CUP 60 ml - graduated



Product Code	26620
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2830724
CND	Y150915
EAN/UPC	8023279266207
GTIN	08023279266207
GMDN	42893

Description

Medicine cup - 60 ml - graduated

- Dimension: 53 x h 53 mm
- Autoclavable 121°C

Stainless Steel AISI 304 holloware with minimum thickness 0.6 mm