



|               |                |
|---------------|----------------|
| Product Code  | 26621          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |     |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

|            |                |
|------------|----------------|
| RDM (NSIS) | 2740318        |
| CND        | V0402          |
| EAN/UPC    | 8023279266214  |
| GTIN       | 08023279266214 |
| GMDN       | 42893          |

## Description

S/S AISI 304 - 18/8 GALLIPOTS - AUTOCLAVABLE 121°C

- ROUND GALLIPOTS
- Ø 208xh 84 mm
- Capacity 1940 ml

For all wash Basins and Gallipots we indicate upper diameter