

GALLIPOT/LOTION BOWL 80 mm - plastic - graduated 200 ml



| | |
|--------------|-------|
| Product Code | 26631 |
|--------------|-------|

| | |
|--------------|------|
| Unit of sale | 1 pc |
|--------------|------|

| | |
|---------------|----|
| Minimum order | 10 |
|---------------|----|

| | |
|------|-------------------|
| Type | No medical device |
|------|-------------------|

| | |
|-------|--|
| Class | |
|-------|--|

| | |
|--------|--------------|
| UK-REP | Not required |
|--------|--------------|

| | |
|--------|--------------|
| CH-REP | Not required |
|--------|--------------|

| | |
|---------|---------------|
| EAN/UPC | 8023279266313 |
|---------|---------------|

| | |
|-------|-----------|
| PARAF | 903162927 |
|-------|-----------|

Description

PLASTIC HOLLOWARE - Autoclavable hundreds of times

Gallipot/lotion bowl
Size: 80 x h 60 mm
Graduation: 200 ml

For the protection and delivery of sterile instruments to operating rooms. For the administration of procedures, holding cleansing solutions and the administration of medicines and fluids.

4 good reasons to choose plastic line:

- 50% lighter than stainless steel
- autoclavable: 1000 cycles guaranteed
- attractive modern design
- cheaper than equivalent stainless steel

Made of polypropylene material, easy to clean, resistant to most healthcare chemicals. Machine, ultrasonic or manual cleaning are permissible. Autoclave sterilization at 134.6°C (+2°C) for a minimum of 2 minutes. EO Ethylene Oxide or Hydrogen Peroxide sterilization both at max. 65°C.

Made in the UK.