



Product Code	26635
Unit of sale	1 pc
Minimum order	10
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2637609
CND	V0499
EAN/UPC	8023279266351
PARAF	903162954
GMDN	42893

## Description

PLASTIC HOLLOWARE - Autoclavable hundreds of times

Kidney Dish  
 Size mm: 200 x h 45  
 Graduation: 500 ml

For the protection and delivery of sterile instruments to operating rooms. For the administration of procedures, holding cleansing solutions and the administration of medicines and fluids.

4 good reasons to choose plastic line:

- 50% lighter than stainless steel
- autoclavable: 1000 cycles guaranteed
- attractive modern design
- cheaper than equivalent stainless steel

Made of polypropylene material, easy to clean, resistant to most healthcare chemicals. Machine, ultrasonic or manual cleaning are permissible. Autoclave sterilization at 134.6°C (+2°C) for a minimum of 2 minutes. EO Ethylene Oxide or Hydrogen Peroxide sterilization both at max. 65°C.

Made in the UK.