

S/S ECONOMIC KIDNEY DISH 12" - 309x149x59 mm



Product Code	26644
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2739353
CND	V0402
EAN/UPC	8023279266443
GTIN	08023279266443
GMDN	42893

Description

Stainless Steel holloware, autoclavable 121°C, with 0.5 mm thickness.

Economic kidney dish - with curved edges - tray 12"

- Size: 309 x 149 x 59 mm
- Capacity: 1500 ml

Technical Specifications

- Autoclavable 121°C
- AISI 304