

S/S PERFORATED INSTRUM. TRAY WITH LID - 223x126x45 mm



Product Code	26646
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2740255
CND	V0402
EAN/UPC	8023279266467
GTIN	08023279266467
GMDN	12143

Description

Instrument tray with angled edges

- Perforated
- With cover
- Dimension: 223 x 126 x h 45 mm
- Autoclavable 121°C

Stainless Steel AISI 304 holloware with minimum thickness 0.6 mm