

MINI INSTRUMENT BAG - nylon black



| | |
|--------------|-------|
| Product Code | 26649 |
|--------------|-------|

| | |
|--------------|------|
| Unit of sale | 1 pc |
|--------------|------|

| | |
|---------------|---|
| Minimum order | 1 |
|---------------|---|

| | |
|------|-------------------|
| Type | No medical device |
|------|-------------------|

| | |
|-------|--|
| Class | |
|-------|--|

| | |
|--------|--------------|
| UK-REP | Not required |
|--------|--------------|

| | |
|--------|--------------|
| CH-REP | Not required |
|--------|--------------|

| | |
|---------|---------------|
| EAN/UPC | 8023279266498 |
|---------|---------------|

| | |
|-------|-----------|
| PARAF | 903119612 |
|-------|-----------|

Description

Mini Instrument bag, with zip closure.

- Material: made in 600D polyester (external)
- Size: 18 x 14 x h 2 cm
- Colour: Black
- For 8 instrum.