

S/S PERFORATED MAYO TRAY - 350x252x16 mm



Product Code	26651
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2740257
CND	V0402
EAN/UPC	8023279266511
GTIN	08023279266511
GMDN	12143

Description

Mayo tray - shallow

- Perforated
- Dimension: 350 x 252 x h 16 mm
- Autoclavable 121°C

Stainless Steel AISI 304 holloware with minimum thickness 0.6 mm