


ALUMINIUM BOX - 17.5X7.6X2 cm



| | |
|---------------|----------------|
| Product Code | 26662 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | No |
| CH-REP | Yes |

| | |
|------------|---------------|
| RDM (NSIS) | 2740347 |
| CND | V0499 |
| EAN/UPC | 8023279266627 |
| GMDN | 13730 |

Description

ALUMINIUM BOX 17.5x7.6xh 2 cm

Useful containers for the sterilization of surgical instruments and to maintain them sterile.
Silicone tray not included.
Nice aluminium box with blue base and silver cover.