



Product Code	26674
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2805729
CND	S010401
EAN/UPC	8023279266740
PARAF	903163069
GTIN	08023279266740
GMDN	13730

## Description

Round sterilized drums for autoclave and hot air sterilizer.

- Diameter: 150 mm
- Height: 150 mm
- Material: Stainless Steel AISI304

Cylindrical containers used to sterilize dressing materials (gauze compress or cotton etc.) in a steam sterilizer (autoclave) or hot air sterilizer and to keep them as "sterile.

Drums have an effective closing lid with a clip lock system, a carrying handle, air vents system to allow steam to circulate freely during the sterilization cycle/process.

Easy to operate lateral air vent system to be manually closed after sterilization.