



|               |                |
|---------------|----------------|
| Product Code  | 26677          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |     |
|--------|-----|
| UK-REP | Yes |
| CH-REP | No  |

|            |                |
|------------|----------------|
| RDM (NSIS) | 2805786        |
| CND        | S010401        |
| EAN/UPC    | 8023279266771  |
| GTIN       | 08023279266771 |
| GMDN       | 13730          |

## Description

Round sterilized drums for autoclave and hot air sterilizer.

- Diameter: 190 mm
- Height: 240 mm
- Material: Stainless Steel AISI304

Cylindrical containers used to sterilize dressing materials (gauze compress or cotton etc.) in a steam sterilizer (autoclave) or hot air sterilizer and to keep them as "sterile".

Drums have an effective closing lid with a clip lock system, a carrying handle, air vents system to allow steam to circulate freely during the sterilization cycle/process.

Easy to operate lateral air vent system to be manually closed after sterilization.