



| | |
|---------------|----------------|
| Product Code | 26684 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | No |

| | |
|------------|----------------|
| RDM (NSIS) | 2941509 |
| CND | S010401 |
| EAN/UPC | 8023279266849 |
| GTIN | 08023279266849 |
| GMDN | 13730 |

Description

Round sterilized drums for autoclave and hot air sterilizer.

- Diameter: 290 mm
- Height: 240 mm
- Material: Stainless Steel AISI304

Cylindrical containers used to sterilize dressing materials (gauze compress or cotton etc.) in a steam sterilizer (autoclave) or hot air sterilizer and to keep them as "sterile".

Drums have an effective closing lid with a clip lock system, a carrying handle, air vents system to allow steam to circulate freely during the sterilization cycle/process.

Easy to operate lateral air vent system to be manually closed after sterilization.