



Product Code	26686
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2805786
CND	S010401
EAN/UPC	8023279266863
PARAF	903163121
GTIN	08023279266863
GMDN	13730

Description

Round sterilized drums for autoclave and hot air sterilizer.

- Diameter: 390 mm
- Height: 250 mm
- Material: Stainless Steel AISI304

Cylindrical containers used to sterilize dressing materials (gauze compress or cotton etc.) in a steam sterilizer (autoclave) or hot air sterilizer and to keep them as "sterile".

Drums have an effective closing lid with a clip lock system, a carrying handle, air vents system to allow steam to circulate freely during the sterilization cycle/process.

Easy to operate lateral air vent system to be manually closed after sterilization.