

# FORCEPS KELLY - straight - 16 cm



|               |                |
|---------------|----------------|
| Product Code  | 26717          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |     |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

|            |               |
|------------|---------------|
| RDM (NSIS) | 1835323       |
| CND        | L031306       |
| EAN/UPC    | 8023279267174 |
| UMDNS      | 11774         |
| PARAF      | 903119776     |
| GMDN       | 62468         |

## Description

Stainless Steel haemostatic Kelly forceps - straight - 16 cm