

# WEITLANER RETRACTOR - sharp - 20 cm



|               |                |
|---------------|----------------|
| Product Code  | 26865          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |     |
|--------|-----|
| UK-REP | No  |
| CH-REP | Yes |

|            |               |
|------------|---------------|
| RDM (NSIS) | 142823        |
| CND        | L9099         |
| EAN/UPC    | 8023279268652 |
| PARAF      | 903603898     |
| GMDN       | 45182         |

## Description

STAINLESS STEEL RETRACTORS

Weitlaner Retractor - sharp 20 cm - 3x4 prongs