


DRESSING KIT 2 - sterile



| | |
|---------------|----------------|
| Product Code | 26929 |
| Unit of sale | 1 kit |
| Minimum order | 1 |
| Type | Medical device |
| Class | I S |

| | |
|------------|---------------|
| RDM (NSIS) | 2353536 |
| CND | T0203 |
| EAN/UPC | 8059020340018 |
| PARAF | 939234783 |

Description

Kit Content:

Haemostatic flow damp
Ring shaped forceps
Metal forceps with plastic ends
Littauer scissors
Gauze swabs 7.5 x 7.5 cm 8 ply
3 Gauze plugs
Gauze swabs 10 x 20 cm - 16 plies
WNW (TNT)+PE drape 32x50 cm
Waste bag