



Product Code	27003
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2391461
CND	V010302
EAN/UPC	8023279270037
PARAF	903164150

Description

"GIMA" DISPOSABLE BLADES - sterile - dispenser box of 100

- Carbon steel

Surgical blades are available in high quality carbon steel or stainless steel.

Manufactured in accordance with requirements of British Standards BS 2982:1992, ISO 7740:1985, ISO 7153-1: 1991, EN 10088-3: 2005.

Each blade is hallmarked and individually packed in aluminium foil.

Carbon Steel blades are VPI paper protected.

Blade size is clearly identified on the reverse side of the foil peel-pack and on product unit boxes.

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, PL, HU, SK, SE, FI, DK, RO, BG, GR, HR, SI, Arabic.

MDR Certified.