

# 3M SURGICAL CLIPPER STARTER KIT (CLIPPER + CHARGER) - 9667L-E



|               |                |
|---------------|----------------|
| Product Code  | 27054          |
| Unit of sale  | 1 kit          |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2248055       |
| CND        | V9099         |
| EAN/UPC    | 8023279270549 |
| PARAF      | 951090051     |

## Description

### 3M™ SURGICAL CLIPPER STARTER KIT (CLIPPER + CHARGER) - 9667L-E

The 3M surgical clipper with its unique pivoting head (when used with 3M clipper blades 9660 code 27055) is a safe, reliable and easy to use device providing a comfortable clipping experience for patients.

Longlasting battery with a short recharge time, a battery indicator light, an ergonomically designed handle and submersibility. What's the result? A better clipping experience for both you and your patient.

Experience the difference a pivoting head makes: When used with the 9660 blade, the unique pivoting head conforms to body contours, making it easy to provide a comfortable clipping experience for patients. You can efficiently cut through all types of head and body hair, wet or dry, in a single pass.

Ergonomic handle with charger.

Single use clipper blade individually sealed.

Clipper blade for head and heavy body hair.

Two type of clipper blades available:

- 3M™ SURGICAL CLIPPER BLADES - 9660 - pivoting (code 27055)
- 3M™ SURGICAL CLIPPER BLADES - 9690 - non pivoting (code 27056)

## Technical Specifications

### Features of the clipper with PIVOTING blades 9660 (code 27055)

- Pivoting head
- Ergonomically designed handle that's easy to hold
- Easy-to-attach blade
- Cuts through all types of head and body hair, wet or dry, in a single pass
- Away-from-the-skin cutting blade minimizes potential for nicks or cuts, which helps reduce the risk of SSIs (surgical site infections)
- Li-ion battery designed to hold a charger longer and battery indicator light
- Approximately 160 minutes of run time on a four-hour recharge
- Submersible for easy cleaning (meets IP X7 rating of submersible in water no deeper than 1 m for up to 30 min)

### Features of the clipper with NON pivoting blades 9690 (code 27056)

- Ergonomically designed handle that's easy to hold
- Easy-to-attach blade
- Quickly and effectively removes hair in trauma and neuro cases
- Away-from-the-skin cutting blade minimizes potential for nicks or cuts, which helps reduce the risk of SSIs (surgical site infections)
- Li-ion battery designed to hold a charger longer and battery indicator light
- Approximately 160 minutes of run time on a four-hour recharge
- Submersible for easy cleaning (meets IP X7 rating of submersible in water no deeper than 1 m for up to 30 min)