

LARGE "POLUS SKAY" MEDICAL BAG



| | |
|---------------|-------------------|
| Product Code | 27108 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | No medical device |
| Class | |

| | |
|---------|---------------|
| UK-REP | Not required |
| CH-REP | Not required |
| EAN/UPC | 8023279271089 |

Description

Rigid bottom, inside zip fastener folding pocket elasticated retainers for instruments on both inside sides.

External dimensions: 45 x 20 x h 25 cm

Made in Italy.