

"UTILITY" DOCTOR BAG - nylon blue



Product Code	27125
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	
UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279271256

Description

UTILITY DOCTOR BAG - blue - empty

Practical doctor bag for easy transport of instruments and medications, provided with 8 internal movable divisions, 2 removable external pockets and 1 removable mini bag (17x3xh 24 cm), 1 fixed external pocket with zip including 4 elastics retaining and other small compartments.

Size: 36x22xh 25 cm.

Made of water resistant 600D polyester
Removable minibag