

*accessori non inclusi
accessories not included
accessoires non inclus*



Product Code	27128
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279271287

Description

PROFESSIONAL BAG - empty

Versatile bag, designed for general practitioners and any professional. It includes: 3 detachable transparent pouches with velcro, 1 ampoule holder with 23 elastic bands with velcro, 1 large pouch with 4 dividers with velcro, clip paper board, document holder, various pockets. Shoulder strap, handle, hidden straps to carry on back.

Made in 70D nylon (external) 600D polyester (internal).

Weight: 3.3 kg.

Maximum load: 5 kg

External size: 42x35x15 cm.