

SMART TROLLEY BAG - medium - blue



Product Code	27141
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279271416

Description

SMART TROLLEY BAG - 56x32.5x49 cm

Trolley bag made of strong, water resistant 600D polyester with durable cold resistant wheels.
Extremely capacious with internal adjustable compartments which permit personalized organization of the space.
Provided with external pockets, internal divisions, and one removable pouch with transparent window.
Waterproof rubber bottom to prevent water infiltration.