

SMART TROLLEY BAG - medium - red



Product Code	27154
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279271546

Description

EMERGENCY SMART TROLLEY BAG - 56X32.5X49 cm

Extremely capacious with internal adjustable compartments which permit personalized organization of the space.

Provided with external pockets, internal divisions, and one removable pouch with transparent window.

Double yellow reflective strips.

Waterproof rubber bottom to prevent water infiltration.