

EMERGENCY BAG - red



Product Code	27165
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8023279271652
PARAF	977811862

Description

EMERGENCY BAG - polyester - red - 35x45xh 21 cm

Capacious bag, despite the compact size. Large central compartment, 1 front and 3 lateral pockets.
May be carried on shoulder (braces) or hand (2 handles).
Made of strong, water resistant material 600D polyester PVC coated.
Provided with yellow reflective strips and rubber bottom.

2 lateral pockets 22x18 cm and 11x18 cm
Lateral pocket 32x18 cm
Front pocket 35x20 cm