

EMERGENCY BAG PVC coated - red



| | |
|---------------|-------------------|
| Product Code | 27166 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | No medical device |
| Class | |

| | |
|---------|---------------|
| UK-REP | Not required |
| CH-REP | Not required |
| EAN/UPC | 8023279271669 |
| PARAF | 977811874 |

Description

EMERGENCY BAG - PVC - red - 35x45xh 21 cm

Capacious bag, despite the compact size. Large central compartment, 1 front and 3 lateral pockets.
May be carried on shoulder (braces) or hand (2 handles).
Made of strong, water resistant material 600D polyester PVC tarpaulin.
Provided with yellow reflective strips and rubber bottom.

2 lateral pockets 22x18 cm and 11x18 cm
Lateral pocket 32x18 cm
Front pocket 35x20 cm