

E4 RECTANGULAR POUCH long with window - yellow



| | |
|---------------|-------------------|
| Product Code | 27194 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | No medical device |
| Class | |

| | |
|---------|---------------|
| UK-REP | Not required |
| CH-REP | Not required |
| EAN/UPC | 8023279271942 |

Description

Pouche provided with a large strip of female velcro on back.
Made in 600D polyester.

Model: E4
Size: 33 x 15 x 11 cm
Colour: yellow
Standard with: 27151/2/7/8, 27175