

THERMAL BAG - yellow nylon



Product Code	27209
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279272093
---------	---------------

PARAF	982516534
-------	-----------

Description

THERMAL BAG - yellow - nylon

Crafted, durable and insulated thermal bag with detachable shoulder strap.
Supplied with pockets and 3 temperature packs, which maintain cold and hot temperature for long period.

External size: 34x25x24 cm - Internal size: 30x23x21 cm

Made in Italy.