

PROFESSIONAL VACCINE CARRIER



Product Code	27211
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8000303301608

Description

PROFESSIONAL VACCINE CARRIER

Compact, extra strong with special insulation by high density PU foam suitable for transportation of vaccine.

Comes complete with 8 empty ice-packs 400 ml suitable for a cold life of 40 hours at 8°C with external temperature of 43°C without opening.

Vaccine storage capacity: 2.6 l.

Body and lid with gaskets for better insulation.

- External size: 28.5 x 25.5 x h 32 cm
- Internal size: 22 x 18 x h 24 cm

Made in Italy.