



Product Code	27218
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	196129
CND	Z12030305
EAN/UPC	7426848150852

## Description

### ENERGY HEAT REGULATING BAG

Heat regulating bag for the warming of IV liquids. Used as an active source of heat is indicated for the prevention and treatment of hypothermia.

Ideal for emergency medicine, hospitals and field camps.

Its compact design allows for use both hung on the stretcher or in the ambulance as a bottle holder or inside bigger containers.

The infusion liquids remain at a constant temperature and the IV liquids can be transported out of the ambulance for a short period.

The main compartment is heated by an electrically fed resistive pad.

Supplied with front pocket for the electric power supply cable, transparent identity pocket, cable with plug, adjustable shoulder strap and manual GB, IT.

Made in Italy.

## Technical Specifications

Padding material: outside PVC, inside nylon 210

Voltage supplies: 12 V DC, absorption 3,75 A

Size: 24x14x33 cm

Weight: 1,2 kg

Car lighter: DIN ISO 4165

Protection: IP3x

Isolation: Class II