



Product Code	27232
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279272321
PARAF	988784233

Description

QUADRANT SMART WATCH

Smart watches designed to monitor different health and sport parameters.

Quadrant for all occasions, from work to leisure.

11 PARAMETERS: BODY TEMPERATURE, HEART RATE, ECG, BPM, SpO2, RESPIRATION RATE, SLEEP, SPORTS, CALORIES, STEPS, WEATHER.

In addition to the above parameters, both watches include stopwatch, music player and photos.
Wireless connection to SIMPLO App, compatible with Android 8.1 and above, iOS 13 and above.

Supplied with power cable and manual in GB, IT, FR, ES, DE available on www.gimaitaly.com.

Technical Specifications

Accelerometer: 3D sensor
 Display: 1.70" IPS TFT
 Data sync: 4.2 BLE
 Temperature range: -10° ~ 50°
 Battery: Lithium polymer
 Battery capacity: 200 MAH
 Battery life stand-by: 30 days
 normal use: 15 days
 frequent use: 7 days
 Charging time: 2-3 hours
 Weight: 55 g
 Water resistance: IP68