



Product Code	27233
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279272338
PARAF	950834224

## Description

### SPORTPIX SMART WATCH

Smart watches designed to monitor different health and sport parameters.

Sportpix is more suitable for sport activities.

11 PARAMETERS: BODY TEMPERATURE, HEART RATE, ECG, BPM, SpO2, RESPIRATION RATE, SLEEP, SPORTS, CALORIES, STEPS, WEATHER.

In addition to the above parameters, both watches include stopwatch, music player and photos.  
Wireless connection to SIMPLO App, compatible with Android 8.1 and above, iOS 13 and above.

Supplied with power cable and manual in GB, IT, FR, ES, DE available on [www.gimaitaly.com](http://www.gimaitaly.com).

## Technical Specifications

Accelerometer: 3D sensor  
 Display: 1.3" OLED  
 Data sync: 4.0 BT  
 Temperature range: -10° ~ 50°  
 Battery: Lithium polymer  
 Battery capacity: 240 MAH  
 Battery life stand-by: 15 days  
 normal use: 5 days  
 frequent use: 3 days  
 Charging time: 2-3 hours  
 Weight: 38 g  
 Water resistance: IP67