



Product Code	27240
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279272406

Description

GIMA DIGITAL BIG DIAL SCALE

Digital scale with extra large, easy to read LCD display (65 mm).
Strong and durable steel construction, comfortable platform.

Lb/kg/st switchable by user.
Supplied in multilingual gift box.

Multilingual manual and box: GB, FR, IT, ES, DE, PT, GR, Arabo

Technical Specifications

Capacity: 225 kg
Sensitivity: 100 g
Minimum measuring weight: 2 kg
Power: lithium battery (included)
Dimension: 430x380x70 mm