

SECA 807 DIGITAL SCALE



Product Code	27257
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	4012030006750

Description

SECA 807 DIGITAL SCALE

Digital scale with advanced technology for enhanced comfort.
Extra-flat design to easily step onto the platform.
Four specially constructed load cells made of stainless steel, guarantee high precision.
Supplied with power sensors, extralarge LCD numbers and high-quality glass surface.
Safety shock resistant glass construction.
Step-off, autohold and automatic switch ON/OFF functions.

Multilingual manual: GB, FR, IT, ES, PT, DE, PL, NL, GR, DK, SE, FI, NO, TR, JP, RU.

Digit height: 29 mm

Technical Specifications

Capacity: 150 kg
Sensitivity: 100 g
Size: 310x310xh 37 mm
Weight: 2 kg
Batteries: 4xAAA 1,5V (not included)