



Product Code	27261
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	599730756514

Description

DIGITAL SMALL PET SCALE

Tray 530x290 mm
 Lb/kg/st switchable by user
 LCD display 25 mm
 Automatic switch-off
 Overload and low battery indication
 Tare and automatic hold functions

Multilingual manual: GB, FR, IT, ES, SE, FI, CZ, RO, GR, HR.

Technical Specifications

Capacity/Sensitivity: 20 kg / 5 g
 Minimum measurable weight: 20 g
 Tray size: 530x290 mm
 Lithium battery: CR2032, 3 V (not included)