



Product Code	27264
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	4006484500605

## Description

### SOEHNLE BABY HEIGHT ROD

Digital and mechanical rod 35 to 80 cm, division 1 cm.  
LCD display, digit height 16 mm.

Manual in 23 languages.

Switchable to inches  
Size: 630x290x70 mm  
Power supply: 2x1.5V AAA batteries

Includes adapter for attachment to baby scales 27306, 27260