



Product Code	27265
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	4012030020770

Description

SECA 813 DIGITAL FLAT SCALE

Electronic scale with high capacity up to 200 kg.
 Extra-wide platform and flat shape facilitate weighing of heavyweight patients and make stepping on very easy.
 Made of solid steel frame with tough rubber coating and non-slip rubber base, ensuring safety and stability.
 Kg/lbs/sts switch.
 Auto-hold, step-off and automatic switch ON/OFF functions.
 Multilingual manual: GB, FR, IT, ES, PT, DE, PL, NL, GR, DK, SE, FI, NO, TR, JP, RU.

Technical Specifications

Capacity: 200 kg - 440 lbs
 Sensitivity: 100 g - 0,2 lbs
 Size: 433x373xh 47 mm
 Platform size: 400x300xh 45 mm
 Weight: 2,9 kg
 Batteries: 4xAAA 1,5V (not included)