



Product Code	27273
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	4006484780915

Description

CHAIR SCALE SOEHNLE 7802

Fully electronic scale: no wear and tear parts, no maintenance.
 Easy to read display with 20 mm high figures.
 Zero setting, taring and hold functions.
 Tiltable arms for easy access of bed-ridden and frail patients.
 Powered by either mains or rechargeable batteries.
 Cable for recharging batteries included.

2 rear swivel castors.

Manual in: DE, GB, FR, NL, ES, IT, PT, CZ, FI, SE, NO, DK, GR, PL, TR.

Made in Germany.

Technical Specifications

Capacity: 200 kg (440 lb)
 Sensitivity: 100 g
 Total weight: 25 kg
 Dimensions: 580x1,100x930 mm
 Power supply: 230 V / 50-60 Hz, 10 mA, 2.2 VA
 Rechargeable batteries