



Product Code	27275
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	4006484781059
PARAF	987476456

## Description

### SOEHNLE 7830

Professional personal scales with stand and gently curved glass indicator. While stepping on the scale, the display is activated automatically and the weight is shown with red backlighted numbers. It does not require any operation and is suitable for medical institutions, wellness areas and pharmacies.

#### Features:

- easy to use, tare key on the back of the display
- glass scale indicator (glass thickness 8 mm)
- bright LED display: 5 digits, height 15 mm

Multilingual manual in: DE, GB, FR, ES, IT, NL, DK, FI, SE, NO, GR, TR, CZ, SK, PL, RU.

Made in Germany.

## Technical Specifications

Capacity: 250 kg  
Sensitivity: 100 g  
Dimension mm: 365x490x890  
Weight: 14.1 kg  
Power Supply: 230 V - 50-60 Hz (available on request 110 V)