



Product Code	27283
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No
RDM (NSIS)	76914
CND	Z12080402
EAN/UPC	4012030008693
GMDN	173918

Description

SECA 384 - electronic

Two scales in one: with a comfortable tray to weigh in a reclining position, easily detachable to be converted into a floor scale. Breast-milk-intake function determines how much milk the baby has suckled.

- Metrological Approval for use in medical field - class III
- Calibration according to class III sensitivity required by 90/384/EEC

Other functions: tare, auto-hold, kg/lbs switch over automatic switch off.
User manual: GB, FR, IT, ES, DE, PT, PL, NL, GR, DK, SE, FL, NO.

Technical Specifications

- Capacity: 20 kg
- Size: 552 x 156 x h 332 mm
- Battery operated
- Sensitivity: 10 g < 10 kg > 20 g
- Weight: 2.8 kg
- Made in Asia