

SECA 756 MECHANICAL SCALE - with height meter - class III



Product Code	27284
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No

RDM (NSIS)	92665
CND	V0399
EAN/UPC	4012030004312
PARAF	987475860
GMDN	173914

Description

SECA 756 - mechanical with height rod (60-200 cm) - class III

Tradition combined with modern design. With a very sensitive precision mechanism to show weight through the aid of a well visible red pointer on the circular dial.

User manual in GB, FR, IT, ES, DE, PT, PL, NL, GR, DK, SE, FI, NO.

Technical Specifications

Capacity: 160 kg
Sensitivity: 500 g
Weight: 18.5 kg
Dimensions: 628x360x1,000 mm
Platform: 360x350 mm with wheels