



Product Code	27291
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No
RDM (NSIS)	2701525
CND	V0399
EAN/UPC	4015672110908
PARAF	985711011

## Description

### OMRON VIVA HBF-222T-EBK

OMRON VIVA is a Medical Device with Clinically Validated Accuracy of key parameters that allow to better understand heart's health and predict cardiovascular events, tracking blood pressure and heart rate.

It monitors 6 precise body metrics: visceral fat, fat, skeletal muscle, resting metabolism, BMI and weight, using Bio-electric Impedance (BIA) Analysis Technology to assess body composition.

Automatic User Recognition function with 4 users + Guest option, 4 seconds fast reading, Auto ON/OFF function and large LCD display.

Results can be shared with the "OMRON connect" App.  
Made of stainless steel and glass.

It includes 4xAAA batteries, setup instructions and multilingual manual: GB, FR, IT, ES, DE, NL, TR, RU, Arabic.

## Technical Specifications

Capacity: 150 kg  
Size: 285x280xh 28 mm  
Weight: 1.6 kg  
Power supply: 4xAAA alkaline batteries (included)  
Memory: 30 readings  
Protection: IP21