

SECA 700 MECHANICAL SCALE - kg/lbs - with height meter



Product Code	27295
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No

RDM (NSIS)	2129637
CND	V0399
EAN/UPC	4012030018838

Description

Mechanical scales with height rod

- SECA 700 - kg/lbs

Mechanical column scale with weight adjustment at eye level.

Weight can be measured in a conventional way: by aligning two chrome plated indicators.

User manual in GB, FR, IT, ES, DE, PT, PL, NL, GR, DK, SE, FI, NO.

Technical Specifications

Capacity: 160 kg - 400 lbs
Sensitivity: 100 g - 1/4 lbs
Dimensions: 520 x 1,556 x 520 mm
Platform: 335 x 80 x 345 mm
Weight: 16.5 kg