

SECA 711 MECHANICAL SCALE - class III - with height meter



Product Code	27296
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No

RDM (NSIS)	92665
CND	V0399
EAN/UPC	4012030003483
PARAF	900425036
GMDN	173936

Description

SECA 711 - class III

Mechanical column scale with weight adjustment at eye level.

Weight can be measured in a conventional way: by aligning two chrome plated indicators.

User manual in GB, FR, IT, ES, DE, PT, PL, NL, GR, DK, SE, FI, NO.

Technical Specifications

Capacity: 220 kg

Sensitivity: 100 g

Dimensions: 520 x 1,556 x 520 mm

Platform: 335 x 80 x 345 mm

Weight: 16.5 kg