

SECA 700 MECHANICAL SCALE - kg - with height meter



Product Code	27297
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No

RDM (NSIS)	1404824
CND	Y0399
EAN/UPC	4012030018791
GMDN	173938

Description

SECA 700 - kg

Mechanical column scale with weight adjustment at eye level.
weight can be measured in a conventional way: by aligning two chrome plated indicators.

User manual in GB, FR, IT, ES, DE, PT, PL, NL, GR, DK, SE, FI, NO.

Technical Specifications

Capacity: 220 kg
Sensitivity: 50 g
Dimensions: 520 x 1,556 x 520 mm
Platform: 335 x 80 x 345 mm
Weight: 16.5 kg