



|               |                |
|---------------|----------------|
| Product Code  | 27299          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I M            |

|        |    |
|--------|----|
| UK-REP | No |
| CH-REP | No |

|            |               |
|------------|---------------|
| RDM (NSIS) | 1403040       |
| CND        | Y0399         |
| EAN/UPC    | 4012030013567 |

## Description

The sturdy classic: with its softly curved weighing tray, the stable all-steel casing and the precise, smooth running weights, this model is the indestructible original baby scale.

User manual: GB, FR, IT, ES, OE, PT, PL.

## Technical Specifications

- Capacity: 16 kg
- Sensitivity: 5 g
- Size: 550 x 180 x h 290 mm
- Weight: 6.4 kg
- Made in Asia