



|               |                   |
|---------------|-------------------|
| Product Code  | 27300             |
| Unit of sale  | 1 pc              |
| Minimum order | 1                 |
| Type          | No medical device |
| Class         |                   |

|         |               |
|---------|---------------|
| UK-REP  | Not required  |
| CH-REP  | Not required  |
| EAN/UPC | 8023279273007 |

Description

GIMA Mechanical baby scale

High quality sturdy mechanical baby scale for accurate and sensitive weighing.  
Supplied with curved weighing tray, aluminium levers and pins, stable casing for safe weighing.  
User manual in GB, FR, IT, ES, PT, DE, RO, GR.

Technical Specifications

Capacity: 20 kg  
Sensitivity: 10 g  
Dimensions: 570x330x285 mm  
Weight: 6.5 kg